

Parametre:

Duette^{HD}

CONTACT LENSES

Material	petrafocon A (center) hemlarafilcon A (skirt)
Dk	Center: 130 Skirt: 84
Water Content	32% (silicone hydrogel skirt)
Diameter	14.5 mm
Base Curves	6.9 to 8.3 in 0.10 mm steps
Skirt Curves	Steep, Medium, Flat, Flat2
Sphere Powers	+20.00D to +8.50D in 0.50D steps +8.00 to -8.00D in 0.25D steps -8.50 to -20.00D in 0.50D steps
Wear indications	Daily wear
Recommended Replacement	6 months

Duette^e

PROGRESSIVE
CONTACT LENSES

Material	petrafocon A (center) hemlarafilcon A (skirt)
Dk	Center: 130 Skirt: 84
Water Content	32% (silicone hydrogel skirt)
Diameter	14.5 mm
Base Curves	7.1 to 8.3 in 0.10 mm steps
Skirt Curves	Medium, Flat, Flat2
Sphere Powers	+5.50D to -8.00D in 0.25D steps -8.50 to -10.00D in 0.50D steps
Add zones	+1,00D, +1,75D, +2,50D
Wear indications	Daily wear
Recommended Replacement	6 months



Ormhojgårdvej 16 · DK-8700 Horsens
Tlf. +45 75 68 10 90 · Fax +45 75 68 10 81
con-lens@con-lens.dk · www.con-lens.dk