

Synergeyes A

Base Curves: 7.00 to 8.20 in 0.10mm steps
Skirt Curves: Base curve +1.00mm or +1.30mm
Sphere Powers: -8.00 to +8.00 in 0.25D steps / +8.50 to +20.00 and -8.50 to -20.00 in 0.50D steps

Synergeyes A Multifocale

Base Curves: 7.00 to 8.20 in 0.10mm steps
Skirt Curves: Base curve +1.00mm or +1.30mm
Sphere Powers: +5.00 to -8.00 in 0.25D steps/ -8.50 to -20.00 in 0.50D steps
Add Powers: +0.75. +1.25. +1.75. +2.25 Segment size: 1.90mm. 2.20mm

ClearKone

Vault Values: 50 μ to 700 μ in 50 μ steps
Skirt Curves: Steep2/Steep/Medium/flat
Sphere Powers: +2.00 to -8.00 in 0.25D steps/ -8.50 to -20.00 and +2.50 to +10.00 in 0.50D steps

Synergeyes KC

Base Curves: 5.30 to 7.10 in 0.20mm steps
Skirt Curves: 5.30 to 6.50, skirt 7.90/8.20/8.50
6.70 and 6.90, skirt 8.20/8.50/8.80
7.10 skirt 8.50/8.80/9.10
Sphere Powers: +8,00 to -20.00D in 0.50 steps

Synergeyes PS

Base Curves: 7.20 to 9.00 in 0.20 steps
Skirt Curves: 8.30 and 8.60. Lift 1. 2. or 3
Sphere Powers: +6.00 to -8.00 in 0.25D steps / -8,5 to -20,00 in 0,50D steps

Diameter: 14.50mm

Materials: Dk 100 Paflucocon D. RGP center / 27% hem-iberfilcon A. soft skirt

Class II UV blocker : >80% UVA and 95% UVB radiation

Daily wear: Replacement at 6 months

