

## Ultrahealth FC parametre:

### Lens Materials:

- 130 Dk rigid gas permeable center
- 84 DK silicone hydrogel soft skirt
- Class II UV blocker > 80% of UVA and >95% of UVB radiation

### Diameter:

14.5mm

### Skirt Curves:

Flat (8.4), Medium (8.1), Steep (7.9)

### Lens Powers:

+10.00 to +2.50 in 0.50D Steps

+2.00 to -8.00 in 0.25D Steps

-8.50 to -20.00 in -0.50D Steps

### Vault Values:

55, 105, 155, 205, 255, 305, 355, 405, 455, 505, 555

### Wear & Replacement:

Daily Wear. Replace at 6 months.

