

Ultrahealth parametre:

Lens Materials:

- 130 Dk rigid gas permeable center
- 84 DK silicone hydrogel soft skirt
- Class II UV blocker > 80% of UVA and >95% of UVB radiation

Diameter:

14.5mm

Skirt Curves:

Flat2, Flat, Medium, Steep

Flat2 only 50-250 vault

Lens Powers:

+10.00 to +2.50 in 0.50D Steps

+2.00 to -8.00 in 0.25D Steps

-8.50 to -20.00 in -0.50D Steps

Vault Values:

50, 100, 150, 200, 250, 300, 350, 400, 450, 500, 550

Wear & Replacement:

Daily Wear. Replace at 6 months.

